

Product Datasheet

(Mouse) Dppa3 Antibody (C-term) (orb1167293)

| | |
|----------------------------|---|
| Description | (Mouse) Dppa3 Antibody (C-term) |
| Species/Host | Rabbit |
| Reactivity | Mouse |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Expiration Date | 12 months from date of receipt. |