

Product Datasheet

Bi-Phospho-MET/HGFR(Y1234/Y1235) Antibody (orb1168225)

| | |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Description | Bi-Phospho-MET/HGFR(Y1234/Y1235) Antibody |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Monoclonal |
| Expiration Date | 12 months from date of receipt. |