

Product Datasheet

UNLB Goat F(ab')₂ anti-Rabbit IgG(H+L), mouse/rat/human ads (orb1168905)

Description	UNLB Goat F(ab') ₂ anti-Rabbit IgG(H+L), mouse/rat/human ads
Conjugation	Unconjugated
Tested Applications	FACS
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.