

Product Datasheet

UNLB Swine anti-Rabbit IgG (H+L), human/rat/mouse/Goat SP ads (orb1168916)

Description	UNLB Swine anti-Rabbit IgG (H+L), human/rat/mouse/Goat SP ads
Conjugation	Unconjugated
Tested Applications	FACS
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.