

Product Datasheet

UNLB Donkey anti-Rabbit IgG(H+L), human/rat/mouse SP ads (orb1168919)

| | |
|----------------------------|---|
| Description | UNLB Donkey anti-Rabbit IgG(H+L), human/rat/mouse SP ads |
| Conjugation | Unconjugated |
| Tested Applications | FACS |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |