

## Product Datasheet

### mouse CD19 Antibody (orb1172733)

<b>Description</b>	Rat monoclonal antibody to mouse CD19
<b>Species/Host</b>	Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	FC
<b>Target</b>	CD19
<b>Preservatives</b>	0.01M TBS (pH7.4)
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rat IgG2a,k
<b>Clonality</b>	Monoclonal
<b>MW</b>	90 kDa
<b>Uniprot ID</b>	<a href="#">P15391</a>
<b>Dilution Range</b>	Flow-Cyt=1:50-100
<b>Expiration Date</b>	12 months from date of receipt.