


Product Datasheet

CGA/HCG- α antibody (orb1172787)

Description	CGA/HCG- α antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Target	CGA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	13 kDa
Uniprot ID	P01215
Dilution Range	WB=1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1:PTM-5848+ competitiv...