

Product Datasheet

Cy3, Goat Anti-Rabbit IgG (orb1173260)

Description	Cy3 conjugated Goat Anti-Rabbit IgG (H+L) secondary antibody
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF
Immunogen	Rabbit IgG whole molecule
Preservatives	Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative, 1% rAlbumin as stablizer and 50% Glycerol.
Form/Appearance	Liquid solution
Storage	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.
Note	For research use only
Application notes	Optimal working dilutions should be determined experimentally by the investigator. Suggested 1:50-1:1000 for most fluorescent applications.
Isotype	Goat IgG
Clonality	Polyclonal
Hazard Information	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)