
Product Datasheet

HRP, Goat Anti-Rabbit IgG (orb1173275)

Description

HRP conjugated Goat Anti-Rabbit IgG (H+L)

Species/Host	Goat
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, IHC-P, WB
Immunogen	Rabbit IgG whole molecule
Preservatives	Liquid in PBS, pH 7.4, containing 1% rAlbumin as stabilizer and 50% Glycerol.
Form/Appearance	Liquid solution
Storage	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.
Note	For research use only
Application notes	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:10000-1:100000), IHC-P (1:500-1:2000), ELISA (1:5000-1:100000).
Isotype	Goat IgG
Clonality	Polyclonal
Hazard Information	The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.
Expiration Date	12 months from date of receipt.