
Product Datasheet

HA Tag Mouse Monoclonal Antibody (4F6) (orb1173282)

Description
Featured Anti-HA Tag Mouse

designed for your tagged protein quan...

Species/Host

Mouse

Conjugation

Unconjugated

Tested Applications

IF, IP, WB

Immunogen

Synthetic Peptide

Preservatives

Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Form/Appearance

Liquid solution

Concentration

1 mg/ml

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:10000), IF (1:2000), IP (1:400).

Isotype

Mouse IgG

Clonality

Monoclonal

Hazard Information

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

Expiration Date

12 months from date of receipt.