
Product Datasheet

GFP Tag Rabbit Polyclonal Antibody (orb1173285)

Description

Featured Anti-GFP Tag Rabbit Polyclonal

tagged protein quantifi...

Species/Host

Rabbit

Conjugation

Unconjugated

Tested Applications

IF, WB

Immunogen

Recombinant Protein

Preservatives

Liquid in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Form/Appearance

Liquid solution

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only

Application notes

Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:5000); IF(1:500).

Isotype

Rabbit IgG

Clonality

Polyclonal

Hazard Information

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

Expiration Date

12 months from date of receipt.