
Product Datasheet

**Magnetic Beads Conjugated Anti-DDDDK Tag
Mouse Monoclonal Antibody (1B10) (orb1173290)**

Description

Featured Anti-DDDDK Tag Mouse

Beads, specially designed for immunoprec...

Species/Host

Mouse

Conjugation

Magnetic Beads

Tested Applications

IP

Immunogen

Synthetic Peptide

Preservatives

Suspended in PBS, pH 7.4, containing 0.02% Sodium Azide as preservative.

Form/Appearance

Liquid solution

Storage

Store at 4°C for 12 months. Avoid freeze-thaw or centrifugation.

Note

For research use only

Application notes

Optimal working dilutions should be determined experimentally by the investigator.

Isotype

Mouse IgG1

Clonality

Monoclonal

Hazard Information

The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.

Expiration Date

12 months from date of receipt.