

Product Datasheet

Human IFI16 ELISA Kit (orb1173603)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 0.16-10 ng/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | Q16666 |
| Sensitivity | 0.062 ng/mL |
| Expiration Date | Please enquire. |

