

Product Datasheet

Human MS4A4A ELISA Kit (orb1173605)

| | |
|----------------------------|---|
| Description | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| Reactivity | Human |
| Range | 15.63-1000 pg/mL |
| Tested Applications | ELISA |
| Note | For research use only |
| Sample Types | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids |
| Uniprot ID | Q96JQ5 |
| Sensitivity | 6.4 pg/mL |
| Expiration Date | Please enquire. |

