
Product Datasheet

Human CD93 Protein, mFc Tag (orb1173679)

Description

Recombinant human CD93 Protein with

| | |
|------------------------|---|
| Reactivity | Human |
| Conjugation | Unconjugated |
| Target | CD93 |
| Preservatives | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions. |
| Form/Appearance | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Storage | Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature. |
| Tag | C-Mouse Fc Tag |
| Note | For research use only |
| Purity | The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining. |
| Source | Mammalian |
| MW | The protein has a predicted molecular mass of 84.5 kDa after removal of the signal peptide. The apparent molecular mass of CD93-mFc is approximately 100-130 kDa due to glycosylation. |
| Uniprot ID | Q9NPY3 |
| Expiration Date | 6 months from date of receipt. |