



Product Datasheet

Human TSLP Protein, hFc Tag (orb1173827)

Description	Recombinant Human TSLP with C-terminal human Fc tag
Reactivity	Human
Conjugation	Unconjugated
Target	TSLP
Preservatives	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Form/Appearance	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Storage	Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing). Lyophilized proteins are shipped at ambient temperature.
Tag	Fc tag
Note	For research use only
Purity	The purity of the protein is greater than 95% as determined by SDS-PAGE and Coomassie blue staining.
Source	Mammalian
MW	41.1 kDa
Uniprot ID	Q969D9
Expiration Date	6 months from date of receipt.



Human TSLP Protein, hFc Tag on SDS-PAGE ...



ELISA plate pre-coated by 5 µg/mL (100 ...