

Product Datasheet

CCR8 Antibody (orb1173975)

| | |
|----------------------------|--|
| Description | Anti-CCR8 antibody(DM195), Rabbit mAb |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FC |
| Target | CCR8 |
| Form/Appearance | Liquid |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Application notes | Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution. |
| Isotype | Rabbit IgG |
| Clonality | Recombinant |
| Uniprot ID | P51685 |
| Dilution Range | ELISA 1/5000-10000 Flow Cyt 1/100 |
| Expiration Date | 12 months from date of receipt. |