

Product Datasheet

CLDN18.2 Antibody (orb1173991)

Description	Anti-CLDN18.2 antibody(DM179), Rabbit mAb
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	CLDN18.2
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	Rabbit IgG
Clonality	Recombinant
Uniprot ID	P56856
Dilution Range	ELISA 1/5000-10000 Flow Cyt 1/100
Expiration Date	12 months from date of receipt.