

Product Datasheet

Cynomolgus CD117 / c-kit Protein (orb1176548)

Description	Cynomolgus CD117, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 26 - Thr 52...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD117
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 57.6 kDa. The protein migrates as 70-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_005555329.1
Purity	95%
Source	Cynomolgus CD117, His Tag (orb1176548) is expressed from human 293 cells (HEK293). It contains AA Gln 26 - Thr 521 (Accession # XP_005555329.1).
MW	57.6 kDa
NCBI	XP_005555329.1
Expiration Date	6 months from date of receipt.