
Product Datasheet

Canine HGF R / c-MET Protein (orb1176554)

Description

Canine HGF R,Fc Tag is expressed from

Glu 25 - Leu 935 (Acc...

Reactivity

Canine

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

HGF R

Preservatives

51 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5

Form/Appearance

Powder

Storage

-20°C

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 128.4 kDa. The protein migrates as 35-40 kDa and 110-120 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_001002963.1

Purity

90%

Source

Canine HGF R,Fc Tag (orb1176554) is expressed from human 293 cells (HEK293). It contains AA Glu 25 - Leu 935 (Accession # Q75ZY9-1).

MW

128.4 kDa

Uniprot ID[Q75ZY9](#)**NCBI**[NP_001002963.1](#)**Expiration Date**

6 months from date of receipt.