
Product Datasheet

Biotinylated Human TROP-2 / TACSTD2 (88-274) Protein (orb1176562)

Description

Biotinylated Human TROP-2 (88-

human 293 cells (HEK293). It contains AA ...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Biotin

Target

TROP-2

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His & C-Avi

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus, followed by an Avi tag (). The protein has a calculated MW of 25.1 kDa. The protein migrates as 35-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_002344.2

Purity

90%

Source

Biotinylated Human TROP-2 (88-274), His, (orb1176562) is expressed from human 293 cells (HEK293). It contains AA Thr 88 - Thr 274 (Accession # P09758-1).

MW

25.1 kDa

Uniprot ID**P09758****NCBI****NP_002344.2****Expiration Date**

6 months from date of receipt.