

Product Datasheet

HighFidelity AF555 PCR Labeling Kit (orb1176617)

Description	HighFidelity AF555 PCR Labeling Kit is designed to produce randomly AF555-modified DNA probes by PCR...
Storage	store at -20 °C
Note	For research use only
Application notes	Spectroscopic Propertie: λ_{exc} 555 nm, λ_{em} 572 nm, ϵ 155.0 L mmol ⁻¹ cm ⁻¹ (Tris-HCl pH 7.5).
SMILES	<chem>C(/C=C/c1cn(c(=O)[nH]c1=O)[C@H]1C[C@@H]([C@H](O1)COP(=O)(O)OP(=O)(O)OP(=O)(O)O)NC(=O)CCCCNC(=O)CCCCCN*</chem>
Expiration Date	12 months from date of receipt.