

## Product Datasheet

# Mouse Angiotensin I Converting Enzyme 2 ACE2 ELISA Kit (orb1176641)

<b>Description</b>	Mouse Angiotensin I Converting Enzyme 2 ACE2 ELISA Kit
<b>Range</b>	62.5--4000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	Mouse Angiotensin I Converting Enzyme 2 ACE2
<b>Storage</b>	Store at 2 - 8°C. Keep reconstituted standard and detection Ab at -20 °C
<b>Note</b>	For research use only
<b>Sensitivity</b>	10.0 pg/mL
<b>Expiration Date</b>	Please enquire.