

## Product Datasheet

# Human Complement Factor D CFD ELISA Kit (orb1176655)

<b>Description</b>	Human Complement Factor D CFD ELISA Kit
<b>Range</b>	312 -20,000 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	Human Complement Factor D CFD
<b>Storage</b>	Store at 2 - 8°C. Keep reconstituted standard and detection Ab at -20 °C
<b>Note</b>	For research use only
<b>Sample Types</b>	100 µL final volume, dilution factor varies on samples
<b>Sensitivity</b>	10.0 pg/mL
<b>Expiration Date</b>	Please enquire.