

## Product Datasheet

# Dog/Canine Complement Component 3 C3 ELISA Kit (orb1176673)

<b>Description</b>	Dog/Canine Complement Component 3 C3 ELISA Kit
<b>Range</b>	6.25 -200 ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	Dog/Canine Complement Component 3 C3
<b>Storage</b>	Store at 2 - 8°C. Keep reconstituted standard and detection Ab at -20 °C
<b>Note</b>	For research use only
<b>Sensitivity</b>	0.5 ng/ml
<b>Expiration Date</b>	Please enquire.