

## Product Datasheet

### Canine Haptoglobin (Hpt)ELISA Kit (orb1176727)

<b>Description</b>	Canine Haptoglobin (Hpt)ELISA Kit
<b>Tested Applications</b>	ELISA
<b>Target</b>	Canine Haptoglobin (Hpt)ELISA Kit
<b>Storage</b>	Store at 2 - 8°C. Keep reconstituted standard and detection Ab at -20 °C
<b>Note</b>	For research use only
<b>Sample Types</b>	50 µL final volume, dilution factor varies on samples
<b>Sensitivity</b>	0.3 ng/ml
<b>Expiration Date</b>	Please enquire.