

Product Datasheet

Donkey Anti-Goat IgG (H&L) -HRP (orb1176773)

Description	Goat polyclonal antibody to Donkey Anti-Goat IgG (H&L) -HRP
Species/Host	Goat
Reactivity	Goat
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Full-length Goat IgG (purified)
Target	Donkey Anti-Goat IgG (H&L) -HRP
Form/Appearance	Lyophilized with PBS, pH 7.4, containing 1% rAlbumin and 0.01% thimerosal
Storage	The antibody is stable in lyophilized form if stored at -20°C. After constitution, store at 4°C for 2-3 weeks. Store at -20°C for a longer time. Aliquot to avoid repeated freezing and thawing
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA: 1:10,000-1:100,000. . Western blot: 1:1,000-1:20,000. Note: The optimal working dilution should be determined by the end user. The suitability of this antibody for other applications not listed here has not been determined but should not be excluded
Expiration Date	12 months from date of receipt.