

Product Datasheet

Goat Anti Rabbit IgG (H&L) (orb1176776)

Description	Goat polyclonal antibody to Goat Anti Rabbit IgG (H&L)
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Full-length rabbit IgG (purified)
Target	Goat Anti Rabbit IgG (H&L)
Form/Appearance	1ml in PBS, pH 7.4, containing 0.01% thimerosal
Storage	Store at 4°C for 2-3 weeks. Store at -20°C for a longer time. Aliquot to avoid repeated freezing and thawing
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	ELISA capture. : 1:100-1:500. . ELISA detection:. 1:10,000-1:30,000. . Western blot. : 1:1,000-1:10,000. Note: The optimal working dilution should be determined by the end user. The suitability of this antibody for other applications not listed here has not been determined but should not be excluded
Expiration Date	12 months from date of receipt.