

Product Datasheet

SMN1/2 antibody (orb1177299)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to SMN1/2 |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | IHC-P, WB |
| Target | SMN1/2 |
| Preservatives | Antibody is lyophilized with 5 mg rAlbumin, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ , 0.05 mg Na ₃ . *carrier free antibody available upon request. |
| Form/Appearance | Lyophilized |
| Storage | At -20 °C for 12 months, as supplied. Store reconstituted antibody at 2-8 °C for one month. For long-term storage, aliquot and store at -20 °C. Avoid repeated freezing and thawing. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Dilution Range | Western blot: 0.1-0.5µg/ml. Immunohistochemistry(Paraffin-embedded Section): 0.5-1µg/ml |
| Expiration Date | 12 months from date of receipt. |