

## Product Datasheet

### Pan-Cadherin antibody (orb1177559)

<b>Description</b>	Mouse monoclonal antibody to CDH1; CDH2; CDH3
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Gallus, Human, Mouse, Rabbit, Rat, Snake
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC-P, WB
<b>Target</b>	CDH1; CDH2; CDH3
<b>Preservatives</b>	Mouse ascites fluid, 1.2% sodium acetate, 2 mg rAlbumin, with 0.01 mg NaN <sub>3</sub> as preservative. *carrier free antibody available upon request.
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	At -20 °C for 12 months, as supplied. Store reconstituted antibody at 2-8 °C for one month. For long-term storage, aliquot and store at -20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	Western blot: 1-2µg/ml. Immunohistochemistry(Paraffin-embedded Section): 2-4µg/ml
<b>Expiration Date</b>	12 months from date of receipt.