

Product Datasheet

Cortisol Binding Globulin Antibody (orb1178535)

Description	Rabbit polyclonal antibody to SERPINA6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Target	SERPINA6
Preservatives	Antibody is lyophilized with 5 mg rAlbumin, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ , 0.05 mg Thimerosal and 0.05 mg NaN ₃ . *carrier free antibody available upon request.
Form/Appearance	Lyophilized
Storage	At -20 °C for 12 months, as supplied. Store reconstituted antibody at 2-8 °C for one month. For long-term storage, aliquot and store at -20 °C. Avoid repeated freezing and thawing.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	Western blot: 0.1-0.5µg/ml. Immunohistochemistry(Paraffin-embedded Section): 0.5-1µg/ml. Immunocytochemistry: 0.5-1µg/ml
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)