

## Product Datasheet

### Goat anti-Rabbit IgG-iFluor 594 (orb1178832)

<b>Description</b>	Goat anti-Rabbit IgG-iFluor 594 is conjugated with high-activity fluorescence dye.
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	iFluor 594
<b>Tested Applications</b>	FACS, IF
<b>Immunogen</b>	Purified rabbit IgG
<b>Preservatives</b>	Phosphate Buffered Saline containing 1% rAlbumin, 0.01% thimerosal, pH 7.4.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.5 mg/mL
<b>Storage</b>	Store at 4°C. Protect from light.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Immunogen affinity chromatography.
<b>Dilution Range</b>	IF 1:200-1:1000 FACS 1:200-1:1000
<b>Expiration Date</b>	12 months from date of receipt.