

Product Datasheet

Goat anti-Rabbit IgG-iFluor 405 (orb1178839)

Description	Goat anti-Rabbit IgG-iFluor 405 is conjugated with high-activity fluorescence dye.
Species/Host	Goat
Reactivity	Rabbit
Conjugation	iFluor 405
Tested Applications	FACS, IF
Immunogen	Purified rabbit IgG
Preservatives	Phosphate Buffered Saline containing 1% rAlbumin, 0.01% thimerosal, pH 7.4.
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store at 4°C. Protect from light.
Note	For research use only
Clonality	Polyclonal
Purity	Immunogen affinity chromatography.
Dilution Range	IF 1:200-1:1000 FACS 1:200-1:1000
Expiration Date	12 months from date of receipt.