

Product Datasheet

Goat anti-Rabbit IgG (H+L)-HRP (orb1178844)

Description	Goat anti-Rabbit IgG (H+L)-HRP is conjugated with high-activity HRP. For direct ELISA, less than 1% ...
Species/Host	Goat
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, IHC, WB
Immunogen	Purified rabbit IgG
Preservatives	Phosphate Buffered Saline containing 1% rAlbumin, 0.01% thimerosal, and 50% glycerol, pH 7.4.
Form/Appearance	Liquid
Concentration	Lot-dependent, please refer to product label.
Storage	Store at 4°C for at least 6 months. Store at -20°C for long-term preservation. Avoid repeated freeze/thaw cycles. Protect from light.
Note	For research use only
Clonality	Polyclonal
Purity	Immunogen affinity chromatography.
Expiration Date	12 months from date of receipt.