

## Product Datasheet

### Goat anti-Rabbit IgG (H+L) (orb1178845)

<b>Description</b>	Goat anti-Rabbit IgG(H+L) is purified from goat antiserum by immunoaffinity chromatography. It react...
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Purified rabbit IgG
<b>Preservatives</b>	Phosphate Buffered Saline containing 0.05% sodium azide, pH 7.4.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	5 mg/mL
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Immunogen affinity chromatography.
<b>Dilution Range</b>	WB 1:5000-1:10000
<b>Expiration Date</b>	12 months from date of receipt.