

Product Datasheet

Mouse anti-His tag mAb (orb1178850)

Description	Mouse anti-His tag mAb is specific to 6X His tag, which can be applied to several immunoassays, incl...
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, IF, IP, WB
Immunogen	6X His tag fusion protein
Preservatives	Phosphate Buffered Saline containing 0.01% thimerosal, pH 7.4.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Monoclonal
Purity	> 95% (SDS-PAGE)
Expiration Date	12 months from date of receipt.