

## Product Datasheet

# SARS-CoV-2 Trimeric Spike Protein (Delta B1.617.2 Variant), His Tag, HEK293 (orb1178866)

<b>Description</b>	SARS-CoV-2 Trimeric Spike Protein (Delta B1.617.2 Variant), His Tag, HEK293
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Stored at -80°C. This product should be used and stored under sterile condition. Avoid repeated free-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	Liquid, Phosphate buffered saline, pH 7.4.
<b>Expiration Date</b>	6 months from date of receipt.