

## Product Datasheet

# MERS-CoV Nucleocapsid Protein, His-SUMO tag, HEK293 (orb1178876)

<b>Description</b>	MERS-CoV Nucleocapsid Protein, His-SUMO tag, HEK293
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Lyophilized protein should be stored at -20°C. Proteins are stable for up to one year from date of receipt at -20°C. Upon reconstitution, protein aliquots should be stored at -20°C or -80°C. Please avoid repeated freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Application notes</b>	It is recommended to reconstitute the lyophilized protein in sterile H <sub>2</sub> O to a concentration not less than 100 µg/mL and incubate the stock solution for at least 20 min to ensure sufficient re-dissolved.
<b>Expiration Date</b>	6 months from date of receipt.