

## Product Datasheet

# Human anti-SARS-CoV & CoV-2 NP Antibody (IgG), clone 104, Ultra Sensitive (orb1178879)

<b>Description</b>	Human anti-SARS-CoV & CoV-2 NP antibody recognize human SARS Coronavirus Nucleocapsid protein (NP). ...
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Preservatives</b>	Phosphate Buffered Saline containing 0.2% Proclin 300, pH 7.4.
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	1 mg/mL
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95% (SDS-PAGE)
<b>Dilution Range</b>	ELISA 1:5000-20000
<b>Expiration Date</b>	12 months from date of receipt.