

Product Datasheet

Human anti-SARS-CoV & CoV-2 Spike antibody (IgA), clone CR3022 (orb1178882)

Description	Human anti-SARS-CoV-2 Spike Antibody [CR3022] recognize human SARS-CoV and CoV-2 Spike protein with ...
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Preservatives	Phosphate Buffered Saline, pH 7.4.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgA
Clonality	Monoclonal
Purity	> 95% (SDS-PAGE)
Dilution Range	ELISA 1:5000-20000
Expiration Date	12 months from date of receipt.