

Product Datasheet

Human anti-SARS-CoV & CoV-2 NP Antibody (IgA) (orb1178887)

Description	Human anti-SARS-CoV & CoV-2 NP antibody recognizes human SARS Coronavirus Nucleocapsid protein (NP)....
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, IF, WB
Preservatives	Phosphate Buffered Saline containing 0.01% thimerosal, pH 7.4.
Form/Appearance	Liquid
Concentration	Lot dependent
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgA
Clonality	Monoclonal
Dilution Range	ELISA 1:5000-20000, WB/IFA/FACS
Expiration Date	12 months from date of receipt.