

Product Datasheet

Mouse anti-SARS-CoV-2 Spike mAb, clone LGSV201 (orb1178888)

Description	This monoclonal antibody is derived from SARS-CoV-2 VLP immune mouse. Coronavirus Spike protein cond...
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Preservatives	Phosphate Buffered Saline containing 0.01% thimerosal, pH 7.4.
Form/Appearance	Liquid
Concentration	Lot dependent
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1 kappa
Clonality	Monoclonal
Purity	> 98% (SDS-PAGE)
Dilution Range	ELISA 1:5000-20000, WB 1:1000-5000, IFA 1:500-1000, FCAS
Expiration Date	12 months from date of receipt.