

Product Datasheet

Human anti-SARS-CoV & CoV-2 NP Antibody (IgG) (orb1178890)

Description	Human anti-SARS-CoV & CoV-2 NP antibody recognize human SARS Coronavirus Nucleocapsid protein (NP). ...
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Preservatives	Phosphate Buffered Saline containing 0.03% Proclin 300, pH 7.4.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	> 98% (SDS-PAGE)
Dilution Range	ELISA 1:5000-20000,WB 1:1000-5000,IFA 1:500-1000,FACS
Expiration Date	12 months from date of receipt.