


Product Datasheet

Recombinant GDNF (Glial-derived neurotrophic factor), Human, AF (orb1179093)

Description	Recombinant GDNF (Glial-derived neurotrophic factor), Human, AF
Reactivity	Human
Endotoxins	< 0.1 EU per 1 µg of the protein by the LAL method.
Tested Applications	ELISA
Preservatives	The protein was lyophilized from a solution containing 20 mM sodium citrate, 0.2 M NaCl, pH 3.5.
Form/Appearance	Lyophilized
Storage	Lyophilized protein should be stored at -20°C. Upon reconstitution, protein aliquots should be stored at -20°C or -80°C.
Tag	His-tag at the N-terminus
Note	For research use only
Protein Sequence	MSPDKQMAVLPRRERNRQAAAANPENSRGKGRRGQRGKNRGCVLTAIHLNVT DLGLGYETKEELIFRYCSGSCDAAETTYDKILKNLSRNRRLVSDKVGQACCRPIAF DDDLSFLDDNLVYHILRKHSKRRCGCI with polyhistidine tag at the C-terminus.
Purity	> 95% as determined by SDS-PAGE. Ni-NTA chromatography.
Source	Escherichia coli
Entrez	2668
Expiration Date	6 months from date of receipt.



SDS- PAGE analysis of recombinant human ...