

Product Datasheet

Recombinant FGF-8a (Fibroblast growth factor-8a), Human, AF (orb1179110)

Description	Recombinant FGF-8a (Fibroblast growth factor-8a), Human, AF
Reactivity	Human
Endotoxins	< 0.1 EU per 1 µg of the protein by the LAL method.
Tested Applications	ELISA
Preservatives	The protein was lyophilized from a solution containing 1X PBS, pH 8.0.
Form/Appearance	Lyophilized
Storage	Lyophilized protein should be stored at -20°C. Upon reconstitution, protein aliquots should be stored at -20°C or -80°C.
Tag	His-tag at the N-terminus
Note	For research use only
Protein Sequence	MQHVREQLVTDQLSRRLIRTYQLYSRTSGKHVQVLANKRINAMAEDGDPFAKL IVETDTFGSRVVRGAETGLYICMNKKGKLIAKSNGKGKDCVFTEIVLENNYTALQ NAKYEGWYMAFTRKGRPRKGSKTRQHQREVHFMKRLPRGHHTTEQSLRFELN YPPFTRSLRGSQRTWAPEPR with polyhistidine tag at the C-terminus.
Purity	> 98% as determined by SDS-PAGE. Ni-NTA chromatography.
Source	Escherichia coli
Entrez	2253
Expiration Date	6 months from date of receipt.