

## Product Datasheet

# Recombinant FGF-16 (Fibroblast growth factor-16), Human, AF (orb1179123)

<b>Description</b>	Recombinant FGF-16 (Fibroblast growth factor-16), Human, AF
<b>Reactivity</b>	Human
<b>Endotoxins</b>	< 0.1 EU per 1 µg of the protein by the LAL method.
<b>Tested Applications</b>	ELISA
<b>Preservatives</b>	The protein was lyophilized from a solution containing 1X PBS, pH 7.4.
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Lyophilized protein should be stored at -20°C. Upon reconstitution, protein aliquots should be stored at -20°C or -80°C.
<b>Tag</b>	His-tag at the N-terminus
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MAEVGGVFASLDWDLHGFSSSLGNVPLADSPGFLNERLGQIEGKLQRGSPTDF AHLKGILRRRQLYCRTGFHLEIFPNGTVHGTRHDHSRFGILEFISLAVGLISIRGVD SGLYLGMNERGELYGSKKLTRECVFREQFEENWYNTYASTLYKHSDSERQYYVA LNKDGSPREGYRTRKHQKFTHFLPRPV DPSKLPMSR DLFHYR with polyhistidine tag at the C-terminus.
<b>Purity</b>	> 95% as determined by SDS-PAGE. Ni-NTA chromatography.
<b>Source</b>	Escherichia coli
<b>Entrez</b>	<b>8823</b>
<b>Expiration Date</b>	6 months from date of receipt.