

## Product Datasheet

### Apcin-A (orb1179327)

<b>Description</b>	Apcin-A
<b>Form/Appearance</b>	Lyophilized
<b>Storage</b>	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Formula</b>	$C_{10}H_{14}Cl_3N_5O_2$
<b>Purity</b>	> 98%
<b>MW</b>	342.61
<b>CAS Number</b>	[1683617-62-0]
<b>SMILES</b>	<chem>O=C(OCCCN)NC(NC1=NC=CC=N1)C(Cl)(Cl)Cl</chem>
<b>Expiration Date</b>	12 months from date of receipt.