

## Product Datasheet

### BI-4020 (orb1179346)

|                        |                                                                                                                                                                                                                                       |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | BI-4020                                                                                                                                                                                                                               |
| <b>Form/Appearance</b> | Lyophilized                                                                                                                                                                                                                           |
| <b>Storage</b>         | Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles. |
| <b>Note</b>            | For research use only                                                                                                                                                                                                                 |
| <b>Formula</b>         | C <sub>30</sub> H <sub>38</sub> N <sub>8</sub> O <sub>2</sub>                                                                                                                                                                         |
| <b>Purity</b>          | > 98%                                                                                                                                                                                                                                 |
| <b>MW</b>              | 542.68                                                                                                                                                                                                                                |
| <b>CAS Number</b>      | N/A                                                                                                                                                                                                                                   |
| <b>SMILES</b>          | <chem>O=C(C1=CC(C)=NC(C2=C(OCCC[C@@H](C)C3)N(C)N=C2)=C1)/N=C(N4)/N3C5=C4C=CC(CN6CCN(C)CC6)=C5</chem>                                                                                                                                  |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                                                                                                                       |