

## Product Datasheet

### Oxotremorine M iodide (orb1179380)

|                        |   |
|------------------------|---|
| <b>Description</b>     | Oxotremorine M iodide   |
| <b>Form/Appearance</b> | Lyophilized   |
| <b>Storage</b>         | Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles. |
| <b>Note</b>            | For research use only   |
| <b>Formula</b>         | $C_{11}H_{19}IN_2O$   |
| <b>Purity</b>          | > 98%   |
| <b>MW</b>              | 322.19  |
| <b>CAS Number</b>      | [713779]  |
| <b>SMILES</b>          | <chem>O=C1N(CC#CC[N+](C)(C)C)CCC1.[I-]</chem>   |
| <b>Expiration Date</b> | 12 months from date of receipt.   |