

Product Datasheet

Aldosterone D8 (orb1179512)

Description	Aldosterone D8
Form/Appearance	Lyophilized
Storage	Store lyophilized at -20°C, keep desiccated. In lyophilized form, the chemical is stable for 36 months. In solution, store at -20°C and use within 3 months to prevent loss of potency. Aliquot to avoid multiple freeze/thaw cycles.
Note	For research use only
Formula	C ₂₁ H ₂₀ D ₈ O ₅
Purity	> 98%
MW	368.49
CAS Number	[1261254-31-2]
SMILES	<chem>OC([2H])([2H])C([C@]1([2H])[C@]23[C@@]([C@@]4([H])[C@]([C@@]5(C(C([2H])([2H])C4)=C([2H])C(C([2H])([2H])C5)=O)C)([H])[C@](OC3O)([H])C2)([H])CC1)=O</chem>
Expiration Date	12 months from date of receipt.